












Medical Temperature Probes

Compatible Brand	LAUNCH Model	DESCRIPTION	CONNECTOR	EMITTER
YSI400 Series 25℃&2.252KΩ	98ME04GA603	Adult Reusable Skin-urface Probe,12.0mm Skin		
	98ME04GA604	Adult Reusable Esophageal/Rectal Probe		
	98ME04GA632	Paediatric Reusable Skin-urface Probe,6.0mm Skin		
	98ME04GA633	Paediatric Reusable Esophageal/Rectal Probe		
	98ME04GA624	Temp Probe Interconnect Calbe		
YSI700 Series 25℃&6/30KΩ	98ME04GA618	Adult Reusable Skin-urface Probe,12.0mm Skin		
	98ME04GA636	Adult Reusable Esophageal/Rectal Probe		
	98ME04GA637	Paediatric Reusable Skin-urface Probe,6.0mm Skin		
	98ME04GA638	Paediatric Reusable Esophageal/Rectal Probe		

















Medical Temperature Probes

Compatible Brand	LAUNCH Model	DESCRIPTION	CONNECTOR	EMITTER
10KΩ Series 25℃&10KΩ	98ME04GA607	Adult Reusable Skin-urface Probe,12.0mm Skin		
	98ME04GA608	Adult Reusable Esophageal/Rectal Probe		
	98ME04GA634	Paediatric Reusable Skin-urface Probe,6.0mm Skin		
	98ME04GA635	Paediatric Reusable Esophageal/Rectal Probe		
50KΩ Series 25℃&50KΩ	98ME04GA641	Adult Reusable Skin-urface Probe,12.0mm Skin		
	98ME04GA642	Adult Reusable Esophageal/Rectal Probe		
	98ME04GA643	Paediatric Reusable Skin-urface Probe,6.0mm Skin		
	98ME04GA644	Paediatric Reusable Esophageal/Rectal Probe		















Medical Temperature Probes

Compatible Brand	LAUNCH Model	DESCRIPTION	CONNECTOR	EMITTER
Philips 飞利浦 25℃&2.252KΩ	98ME04GA619	Adult Reusable Skin-urface Probe,12.0mm Skin		
	98ME04GA620	Adult Reusable Esophageal/Rectal Probe		
	98ME04GA639	Paediatric Reusable Skin-urface Probe,6.0mm Skin		
	98ME04GA640	Paediatric Reusable Esophageal/Rectal Probe		
	98ME04GA625	Temp Probe Interconnect Calbe		
Mindray 迈瑞 25℃&2.252KΩ	98ME04GA641	Adult Reusable Skin-urface Probe,12.0mm Skin		
	98ME04GA642	Adult Reusable Esophageal/Rectal Probe		
	98ME04GA643	Paediatric Reusable Skin-urface Probe,6.0mm Skin		
	98ME04GA644	Paediatric Reusable Esophageal/Rectal Probe		








IBP Cable

Compatible Brand	LATEST P/N	DESCRIPTION	INSTRUMENT CONNECTOR	Yoke type
德恩 Datex	98ME06GA735	UTAH Transducer Adapter Cable	 Round 10J 	
	98ME06GA736	ABBOTT Transducer Adapter Cable		
	98ME06GA737	Edward Transducer Adapter Cable		
	98ME06GA738	BD Transducer Adapter Cable		
	98ME06GA740	Philips Transducer Adapter Cable		
	98ME06GA739	MEDIX Transducer Adapter Cable		
GE-马葵 GE- MARQUETTE	98ME06GA704	UTAH Transducer Adapter Cable	 Square 11P 	
	98ME06GA703	ABBOTT Transducer Adapter Cable		
	98ME06GA713	Edward Transducer Adapter Cable		
	98ME06GA714	BD Transducer Adapter Cable		
	98ME06GA716	Philips Transducer Adapter Cable		
	98ME06GA715	MEDIX Transducer Adapter Cable		

IBP Cable

Compatible Brand	LATEST P/N	DESCRIPTION	INSTRUMENT CONNECTOR	Yoke type
迈瑞 Mindray	98ME06GA718	UTAH Transducer Adapter Cable	 Round 12P	
	98ME06GA717	ABBOTT Transducer Adapter Cable		
	98ME06GA719	Edward Transducer Adapter Cable		
	98ME06GA720	BD Transducer Adapter Cable		
	98ME06GA722	Philips Transducer Adapter Cable		
	98ME06GA721	MEDIX Transducer Adapter Cable		
飞利浦 PHILIPS	98ME06GA728	UTAH Transducer Adapter Cable	 Square 12P	
	98ME06GA729	ABBOTT Transducer Adapter Cable		
	98ME06GA731	Edward Transducer Adapter Cable		
	98ME06GA730	BD Transducer Adapter Cable		
	98ME06GA733	Philips Transducer Adapter Cable		
	98ME06GA732	MEDIX Transducer Adapter Cable		


IBP Cable

Compatible Brand	LATEST P/N	DESCRIPTION	INSTRUMENT CONNECTOR	Yoke type
太空 Spacelabs 迈瑞 Mindray	98ME06GA718	UTAH Transducer Adapter Cable	 Round 6P	
	98ME06GA717	ABBOTT Transducer Adapter Cable		
	98ME06GA719	Edward Transducer Adapter Cable		
	98ME06GA720	BD Transducer Adapter Cable		
	98ME06GA722	Philips Transducer Adapter Cable		
	98ME06GA721	MEDIX Transducer Adapter Cable		

How is the blood pressure monitor working?

- 1.pierce the canula into blood vessel,ventricle,encephalic or muscle.
- 2.the transducer can tranfer the pressure into signal.
- 3.the liquid between canula and transducer is media of pressure signal.
- 4.the monitor system then magnify the signal and display it.

See illustrations :





Non-Invasive Blood Pressure(NIBP) Cuff & Tubes

Reusable Blood Pressure Cuffs

Compatible Brand	LAUNCH Model	DESCRIPTION	ILLUSTRATIONS
<p>PHILIPS Colin Datascope FUKUDA Spacelabs Welch-Allyn Simens Mindray Goldway Criticare</p>	98ME07GB786	Adult Reusable Blood Pressure Cuff	 <p>Single-tube</p> 
	98ME07GB799	Child Reusable Blood Pressure Cuff	
	98ME07GB802	Infant Reusable Blood Pressure Cuff	
	98ME07GB804	Neonate Reusable Blood Pressure Cuff	
	98ME07GA034	Adult Thigh Reusable Blood Pressure Cuff	
<p>GE-Marquette Datex-Ohmeda Critikon New Welch-Allyn</p>	98ME07GB787	Adult Reusable Blood Pressure Cuff	 <p>Two-tube</p> 
	98ME07GB800	Child Reusable Blood Pressure Cuff	
	98ME07GB803	Infant Reusable Blood Pressure Cuff	
	98ME07GB805	Neonate Reusable Blood Pressure Cuff	
	98ME07GA035	Adult Thigh Reusable Blood Pressure Cuff	

Non-Invasive Blood Pressure(NIBP) Cuff & Tubes

Disposable Blood Pressure Cuffs

Compatible Brand	Limb cir.	DESCRIPTION	ILLUSTRATIONS
PHILIPS Colin Datascope FUKUDA Spacelabs Welch-Allyn Simens Mindray Goldway Criticare	3.3cm-5.6cm	1# Neonate Disposable Blood Pressure Cuff	 <p>Single-tube</p>
	4.2cm-7.1cm	2# Neonate Disposable Blood Pressure Cuff	
	5.0cm-10.5cm	3# Neonate Disposable Blood Pressure Cuff	
	6.9cm-11.7cm	4# Neonate Disposable Blood Pressure Cuff	
	8.9cm-15.0cm	5# Neonate Disposable Blood Pressure Cuff	
	12.4cm-16.8cm	Infant Disposable Blood Pressure Cuff	
	20.0cm-27.0cm	Child Disposable Blood Pressure Cuff	
	25.3cm-34.3cm	Adult Disposable Blood Pressure Cuff	
32.1cm-43.4cm	Large Adult Disposable Blood Pressure Cuff		
GE-Marquette Datex-Ohmeda Critikon New Welch-Allyn	3.3cm-5.6cm	1# Neonate Disposable Blood Pressure Cuff	 <p>Two-tube</p>
	4.2cm-7.1cm	2# Neonate Disposable Blood Pressure Cuff	
	5.0cm-10.5cm	3# Neonate Disposable Blood Pressure Cuff	
	6.9cm-11.7cm	4# Neonate Disposable Blood Pressure Cuff	
	8.9cm-15.0cm	5# Neonate Disposable Blood Pressure Cuff	
	12.4cm-16.8cm	Infant Disposable Blood Pressure Cuff	
	20.0cm-27.0cm	Child Disposable Blood Pressure Cuff	
	25.3cm-34.3cm	Adult Disposable Blood Pressure Cuff	
32.1cm-43.4cm	Large Adult Disposable Blood Pressure Cuff		

Air Hose Connectors:



Non-Invasive Blood Pressure(NIBP) Cuff & Tubes

Pressure Cuff Interconnect Tubing

Compatible Brand	FORMER P/N	LATEST P/N	Compatible Models	ILLUSTRATIONS	
飞利浦 PHILIPS	98ME07GC002	DG147086001	Connects all sizes of Philips' reusable and disposable adult and pediatric NIBP cuffs to the monitor. Use with adult and pediatric cuffs only. CANNOT BE USED WITH NEONATAL CUFFS. Use with all Philips models except A1/A3		
	98ME07GA019	DG085085001	NIBP hose connects to 40401X, M45xx, M157x series pressure cuffs. Use with Philips Instrument: A1/A3, M3921A, M3923A, M3926A, M3927A, M3928A, M3929A		
Simens Mindray Spacelabs MEK, NPB	98ME07GA019	DG085085001	Simens, Mindray, MEK, NPB Spacelabs 90496 Module		
科瑞康 Creative	98ME07GC003	DG146085002			
金科威 Goldway		DG152102001			
GE Dinamap, Pro, MPS		DG153103001	DINAMAP, Pro, Compact, MPS Series		
GE marquette	98ME07GC001	DG146085001	Eagle 3000, Dash 2000/3000/4000, Tram 400/450/451, Tram 800. 850/851. Corometrics 120		
		DG154104001			